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#### Resposta:

```
{
  "success": true,
  "response": "Your heart rate variability (HRV) is looking excellent! Your HRV score is 75.0/100",
  "timestamp": "2025-10-21T12:30:45.123456",
  "session_id": "user_123",
  "tokens_used": 245,
  "model": "gpt-4o-mini",
  "error": null
}
```

### Exemplo com curl:

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "Hello! Can you tell me about my current health status?",
  "session_id": "test_user_001",
  "include_biosignal_context": true
}'
```

## 2. Obter Histórico de Conversas

**GET** `/api/v1/chat/history/{session_id}`

Retrieve the conversation history for a specific session.

### Resposta:

```
{
  "session_id": "user_123",
  "history": [
    {
      "role": "user",
      "content": "How is my HRV?"
    },
    {
      "role": "assistant",
      "content": "Your HRV score is 75.0/100..."
    }
  ],
  "message_count": 2
}
```

### Exemplo com curl:

```
curl https://wearable-biosignal-backend.onrender.com/api/v1/chat/history/test_user_001
```

### 3. Limpar Histórico de Conversas

**DELETE** `/api/v1/chat/history/{session_id}`

Delete all conversation history for a session.

**Resposta:**

```
{
  "success": true,
  "message": "Conversation history cleared for session: user_123"
}
```

**Exemplo com curl:**

```
curl -X DELETE https://wearable-biosignal-backend.onrender.com/api/v1/chat/history/test_user_001
```

## Contexto de Biosinais

When `include_biosignal_context` is set to `true`, LIA receives comprehensive real-time data:

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### Clarity™ Layer

- Overall Quality Score (0-1)
- Signal-to-Noise Ratio (dB)
- Quality Assessment (excellent/good/fair/poor)
- Detected Artifacts

## iFRS™ Layer

- Dominant Frequency (Hz)
- Rhythm Classification (normal\_sinus, elevated, low, irregular, athletic)
- HRV Score (0-100)
- HRV Metrics (RMSSD, SDNN, pNN50)
- LF/HF Ratio
- Respiratory Rate

## Timesystems™ Layer

- Pattern Type (stable, increasing, decreasing, oscillating, irregular)
- Circadian Phase (morning, afternoon, evening, night)
- Temporal Consistency Score (0-1)
- Rhythm Score (0-100)
- Circadian Alignment Score (0-1)

## LIA Health Insights

- Detected Health Condition
- Confidence Level (0-1)
- Overall Pontuação de Bem-Estar (0-100)
- Wellness Breakdown (cardiovascular, respiratory, activity, stress)
- Risk Factors
- Positive Indicators
- Recomendações Personalizadas

## Exemplos de Casos de Uso

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```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "How am I doing right now?",
  "session_id": "user_001",
  "include_biosignal_context": true
}'
```

## 2. Entendendo as Métricas

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "What does my LF/HF ratio mean?",
  "session_id": "user_001",
  "include_biosignal_context": true
}'
```

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```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "Is this a good time to exercise based on my circadian phase?",
  "session_id": "user_001",
  "include_biosignal_context": true
}'
```

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```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "What are some tips to improve my sleep quality?",
  "session_id": "user_001",
  "include_biosignal_context": false
}'
```

O servidor exibirá:

- ✓ LIA Core™ Chat Engine initialized with OpenAI
- ✓ Backend ready to accept connections on <https://wearable-biosignal-backend.onrender.com>

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