

LIA Core™ - Exemplos de Comandos cURL

Gerado em: 21/10/2025

LIA Core™ - Exemplos de Comandos cURL

Guia de referência rápida para testar o Módulo Conversacional LIA Core™ usando comandos curl.

Conversa Básica sobre Estado de Saúde

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "Hello! Can you tell me about my current health status?",
  "session_id": "user_001",
  "include_biosignal_context": true
}'
```

Perguntar Sobre Métricas Específicas

Variabilidade da Frequência Cardíaca (VFC)

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "What does my HRV score mean? Is it good or bad?",
  "session_id": "user_001",
  "include_biosignal_context": true
}'
```

Oxigênio no Sangue (SpO2)

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "How is my blood oxygen level? Should I be concerned?",
  "session_id": "user_001",
  "include_biosignal_context": true
}'
```

Ritmo Circadiano

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "What does my circadian phase indicate? Is this a good time for exercise?",
  "session_id": "user_001",
  "include_biosignal_context": true
}'
```

Pontuação de Bem-Estar

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "Explain my overall wellness score. What areas should I focus on improving?",
  "session_id": "user_001",
  "include_biosignal_context": true
}'
```

Perguntas Gerais de Saúde (Sem Dados Atuais)

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "What are the best ways to improve my sleep quality?",
  "session_id": "user_001",
  "include_biosignal_context": false
}'
```

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "What exercises are best for improving cardiovascular health?",
  "session_id": "user_001",
  "include_biosignal_context": false
}'
```

Gerenciamento de Conversas

Obter Histórico de Conversas

```
curl https://wearable-biosignal-backend.onrender.com/api/v1/chat/history/user_001
```

Limpar Histórico de Conversas

```
curl -X DELETE https://wearable-biosignal-backend.onrender.com/api/v1/chat/history/user_001
```

Obter Dados Atuais de Biosinais

```
# Get processed data from all layers  
curl https://wearable-biosignal-backend.onrender.com/api/v1/stream | python3 -m json.tool
```

```
# Get prediction/condition  
curl https://wearable-biosignal-backend.onrender.com/api/v1/predict | python3 -m json.tool
```

```
# Get processing logs  
curl https://wearable-biosignal-backend.onrender.com/api/v1/logs/processing?limit=50 | python3 -
```

Status do Sistema

```
# Check API info  
curl https://wearable-biosignal-backend.onrender.com/ | python3 -m json.tool
```

```
# Health check  
curl https://wearable-biosignal-backend.onrender.com/api/v1/health | python3 -m json.tool
```

Exemplo de Conversa Multi-turno

```
# Message 1
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "Hi LIA, how am I doing today?",
  "session_id": "conversation_demo",
  "include_biosignal_context": true
}'

# Message 2 (follows context from Message 1)
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "Can you explain my HRV in more detail?",
  "session_id": "conversation_demo",
  "include_biosignal_context": true
}'

# Message 3 (follows context from Messages 1 & 2)
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "What should I do to improve it?",
  "session_id": "conversation_demo",
  "include_biosignal_context": true
}'

# View the full conversation
curl https://wearable-biosignal-backend.onrender.com/api/v1/chat/history/conversation_demo | python3
```

Formatar Respostas JSON

Add `| python3 -m json.tool` to any curl command for formatted output:

```
curl -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "How is my heart rate?",
  "session_id": "user_001",
  "include_biosignal_context": true
}' | python3 -m json.tool
```

Extrair Apenas a Mensagem de Resposta

```
curl -s -X POST https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{
  "message": "How is my health?",
  "session_id": "user_001",
  "include_biosignal_context": true
}' | python3 -c "import sys, json; print(json.load(sys.stdin)['response'])"
```

Comandos Úteis de Uma Linha

```
# Quick health check
curl -s https://wearable-biosignal-backend.onrender.com/api/v1/chat \
-H "Content-Type: application/json" \
-d '{"message":"Quick health summary?","session_id":"quick"}' \
| python3 -c "import sys, json; print(json.load(sys.stdin)['response'])"

# Check if server is running
curl -s https://wearable-biosignal-backend.onrender.com/api/v1/health | python3 -c "import sys, json; print(json.load(sys.stdin))"

# Get current heart rate
curl -s https://wearable-biosignal-backend.onrender.com/api/v1/stream | python3 -c "import sys, json; print(json.load(sys.stdin))"

# Get wellness score
curl -s https://wearable-biosignal-backend.onrender.com/api/v1/stream | python3 -c "import sys, json; print(json.load(sys.stdin))"
```

Documentação Interativa da API

Em vez de curl, você também pode testar a API usando a interface Swagger interativa:

Abrir no navegador: **<https://wearable-biosignal-backend.onrender.com/docs>**

Isso fornece uma interface web onde você pode:

- Ver todos os endpoints disponíveis
- Testar requisições diretamente no navegador
- Visualizar esquemas de requisição/resposta
- Experimentar diferentes parâmetros

Nota: Todos os exemplos assumem que o servidor está rodando em <https://wearable-biosignal-backend.onrender.com>.

LIA Core™ - Sistema de Monitoramento de Biosinais